

NEW!

The ExactaMed™ Oral/ Enteral Dispenser System

Multiple use dispensers for administration of oral liquid medication, feeds or flushing directly into the mouth or via a feeding tube.

CE
0543



ExactaMed™
Oral/Enteral
Dispensers with
female luer and
female luer lock tips
NOW AVAILABLE
in 1, 3, 5, 10 and
20 mL sizes



Choose purple. Choose Baxa.

0845 609 0250 choosepurple@baxa.co.uk www.baxa.com/choosepurple

Baxa 

Product Information

ExactaMed™ Oral/Enteral Dispensers

Multiple use dispensers for administration of oral liquid medication either by mouth or via a feeding tube. Validated for multiple use on same named patient.



Oral Tip

Graduations	Size	Single Order Quantity (pcs)	Order Code	NHS Supplychain Order Code	Tip
0.01 mL	1 mL	100	1701	FTA093	Oral
0.1 mL	3 mL	100	1702	FTA094	Oral
0.2 mL	5 mL	100	1705	FTA095	Oral
0.5 mL	10 mL	100	1710	FTA096	Oral
1 mL	20 mL	50	1720	FTA097	Oral
1 mL	60 mL	30	1790	FTA098	Oral



Female Luer (inset) & Female Luer Lock Tips

Graduations	Size	Single Order Quantity (pcs)	Order Code	NHS Supplychain Order Code	Tip
0.01 mL	1 mL	100	1501	FTA176	FL
0.1 mL	3 mL	100	1502	FTA177	FL
0.2 mL	5 mL	100	1505	FTA178	FLL
0.5 mL	10 mL	100	1510	FTA179	FLL
1 mL	20 mL	50	1520	FTA180	FLL
1 mL	60 mL	30	1760	FTA023	FLL



Catheter Tip

Graduations	Size	Single Order Quantity (pcs)	Order Code	NHS Supplychain Order Code	Tip
1 mL	60 mL	30	1750	FTA022	Catheter

Specifications: Clear barrel, Purple plunger, Individually wrapped, Sterile, Latex free, Labelled Oral/Enteral on packaging and barrel.

Availability: Via NHS Supplychain catalogue or direct from Baxa Ltd.

Warning • **DO NOT** autoclave any Baxa Oral/Enteral Dispenser.

- **REPLACE** the dispenser if:
There is any visible damage to barrel or tip
The dose markings are no longer clear
The plunger seal is damaged
It becomes difficult to operate

European Headquarters

Baxa Ltd
Radius Court
Eastern Road
Bracknell
Berkshire
RG12 2UP
UK
Tel: +44 (0) 1344 392 902
Fax: +44 (0) 1344 392 903
e-mail: info@baxa.co.uk
web: www.baxa.com

Austria

Tel: 0800 296 745
Fax: 0800 296 746

Belgium

Tel: 0800 900 94
Fax: 0800 900 87

Denmark

Tel: 808 876 61
Fax: 808 876 60

Finland

Tel: 0800 914 524
Fax: 0800 914 523

France

Tel: 0800 900 189
Fax: 0800 908 436

Germany

Tel: 0800 182 6496
Fax: 0800 182 5856

Luxembourg

Tel: 800 257 65
Fax: 800 257 63

The Netherlands

Tel: 08000 222 530
Fax: 08000 221 990

Switzerland

Tel: 0800 56 17 94
Fax: 0800 56 17 93